

## EASTERN UNIVERSITY, SRI LANKA FIRST YEAR FIRST SEMESTER EXAMINATION IN SCIENCE-2015/2016 (SEPTEMBER' 2018)-REPEAT CH 151-QUANTITATIVE AND QUALITATIVE INORGANIC ANALYSES [OLD SYLLABUS]

**Answer All Questions** 

**Time: Three Hours** 

- 1. You are provided with a mixture  $\underline{A}$  containing two inorganic cations. *Analyze* the mixture  $\underline{A}$  qualitatively and *record* your observations, inferences and conclusion. *Carryout* one confirmatory test for each identified cation.
- A mixture <u>B</u> contains two inorganic anions. *Perform* the following tests and *record* your observations, inferences and conclusion. *Carryout* one confirmatory test for each identified anion.
  - a) Add dil. H<sub>2</sub>SO<sub>4</sub>, warm and *test* for evolved gas
  - b) *Prepare* an aqueous solution of the given sample and *perform* the following experiments.
    - i) Add dil. HNO<sub>3</sub> and AgNO<sub>3</sub>
    - ii) Add dil. HNO<sub>3</sub> and BaCl<sub>2</sub>
    - Add dil H<sub>2</sub>SO<sub>4</sub> and test the evolved gas with filter paper soaked in Lead acetate.
    - iv) Add CaCl<sub>2</sub> and acetic acid to the solution

1

- Acidify with dil. H<sub>2</sub>SO<sub>4</sub> and add freshly prepared FeSO<sub>4</sub> and few drops of con, H<sub>2</sub>SO<sub>4</sub>.
- 3. You are provided with the following solutions.
  - a) An aqueous solution of  $\underline{\mathbf{X}}$  of 0.01 M HCl.
  - b) An aqueous solution of  $\underline{\mathbf{Y}}$  containing 0.5 g mixture of Na<sub>2</sub>CO<sub>3</sub> and NaOH in one litre.

Pipette 10.0 ml solution  $\underline{Y}$  into a titration flask and titrate it against standard . solution  $\underline{X}$  using methyl red as an indicator (take three readings).

Pipette 10.0 ml solution  $\underline{\mathbf{Y}}$  into a titration flask and titrate it against standard solution  $\underline{\mathbf{X}}$  using phenolphthalein as an indicator (take three readings).

- i) Tabulate your readings.
- ii) Calculate the weight of NaOH and Na<sub>2</sub>CO<sub>3</sub> in the mixture.
- iii) Calculate the percentage composition of NaOH and Na<sub>2</sub>CO<sub>3</sub> in the mixture.

2

## End of paper